

Instructions for Adding the CODE STAT site to your LIFEPAK Device

- On your LIFEPAK device, press **On** and simultaneously hold down the **Options** and **Event** buttons to enter setup mode
- Enter the passcode defined for your system (default is **0000**)
- Select **Transmission**
- Confirm **Default Report** is 12-Lead and **Wireless** is turned OFF
- Select **Sites**
- Add a new SITE and Name it **CODE STAT**(Be sure you enter a space between CODE and STAT)
- Select END and Confirm that the **Output Port** to **Direct Connect** and power off the LIFEPAK device

Instructions for Transmitting to the CODE STAT Site

- Transmit **only Cardiac Arrests** to the CODE STAT Site
- You **MUST** change the report to **ALL** before sending to the CODE STAT Site. (Your devices are currently setup with 12-Lead as the Default Report)

Please contact the LIFENET Smart Desk at 1-800-732-3081 or Shelli Baufield at 763-754-1680 for assistance if needed